# **EVOSEP ONE**

A Separation Tool Designed for Clinical Omics

EVUSEP

## EVOSEP

Evosep aims to improve quality of life and patient care by radically innovating protein based clinical diagnostics.

EVA

### Clinical Proteomics must be Sensitive, Fast, and Robust

- More uptime with improved reliability and robustness
- Increased productivity with higher throughput and better duty cycle utilization
- Increased performance with better data quality



"Robust and fast workflows are indispensable for successfully implementing large-cohort clinical plasma proteomics studies",

#### Professor Matthias Mann,

Max-Planck Institute of Biochemistry, Proteomics and Signal Transduction

### **Optimized standard methods**

Throughput (samples/ day)	Cycle Time (min)	Gradient Length (min)	Flow rate (µl/min)	Column dimensions (length/ID/bead size)
300	4,8	3,2	4,0	4cm / 150µm / 3µm
200	7,2	5,6	2,0	4cm / 150µm / 3µm
100	14,4	11,5	1,5	8cm / 100µm / 3µm
60	24,0	21,0	1,0	8cm / 100µm / 3µm
30	48,0	44,0	0,5	15cm / 100µm / 3µm

#### **High instrumentation efficiency**



## Disposable trap columns

- The sample is loaded and desalted offline on a C18 pipette tip, the Evotip
- The autosampler picks up the tip and places it in the injection port

#### **Partial elution**

- Leaves contaminants on the Evotip disposable trap columns
- · Elutes analytically relevant peptides
- Minimizes cross-contamination
- Extends column lifetime

\*Methods subject to change

#### Low maintenance components

All the elution and gradient formation happen at low pressure, ensuring minimal wear and tear

Low pressure pumps

#### Simplified workflows

Integrating elution with liquid chromatography removes sample handling steps as well as reducing injection cycle overheads



### **Robust throughput**

- · 10x reduction in carry-over
- · 1,000 of runs on each analytical column
- · 30,000 sample service interval



The Broad Institute, Cambridge, MA



1500 injections of 1ug HeLa on the same column

#### **Chromatographic performance**



- Intelligent pump control allows
  efficient use of acquisition
  window
- Peptide mix in HeLa background used to calculate peak data
- Overlaid runs show good reproducibility in complex standard

Throughput	FWHM	Peak capacity
Samples/day	[sec]	4σ / FWMM
300	1.5	35 / 72
200	2.2	41 / 97
100	4.0	63 / 129
60	6.8	79 / 161
30	15.2	111 / 216



#### **Reproducible plasma performance**



Retention time stability of selected peptides over 96 runs Pearson correlation matrix comparing the 96 runs

Low carry-over between runs (<0.1%)

Data courtesy: Dr. Philipp Geyer, Max Planck Inst., Martinsried

### Fast and sensitive HCP screening

- · Bottom-up mass spectrometry protocols are established for analyzing host cell proteins (HCPs)
- Detecting low abundance HCPs confidently and with high throughput is now possible with the Evosep One paired with a trapped ion mobility separation TOF instrument



61 HCPs identified in 21 minutes from 25 μg NISTmAb (according to Huang et al., Anal. Chem. 2017, 89, 5436-5444) equivalent to approximately 1 μg load.

- 3 times faster\*
- 2 times more sensitive\*

\*than 1h industry standard

Data courtesy: FN-06 flash note; Bruker Daltonik GmbH, Germany

### Protein identification with short gradients



- Stable protein and peptide ID rates over 96 injections with the 200 samples/day methods
- Examples of identification rates possible with different methods.

All measurements were performed on a Bruker Daltonics timsTOF Pro with 50 ng of HeLa digest and at least 15 replicates.

Data courtesy: LCMS-141 app note, Bruker Daltonik GmbH, Germany

### **Fast generation of DIA libraries**

Spectral library of mammalian cell line proteomes made from 46 fractions of HeLa digests.

Total time	Gradient time	Peptides	Proteins	
18.4h	16.1h	132,850	9,918	

Example ID and quantitation performance with 60 samples/day method



Data courtesy: Prof. Jesper Olsen and Dorte Bekker-Jensen, Novo Nordisk Center for Proteome Research, Copenhagen on a Thermo Scientific Q Exactive HF-X

### Large clinical proteomics sample set





Data courtesy: Dr. Ben Collins, Dr. Evan Williams and Prof. Ruedi Aebersold, ETH Zürich on a Sciex TripleTOF 6600

Enables large cohort studies in a fast and robust manner with

- Stable retention times
- · High data consistency

across hundreds of injections

#### Contact info:

Evosep Thriges Plads 6 DK-5000 Odense C Denmark

sales@evosep.com

# EV**U**SEP

Making clinical proteomics 100 times more robust and 10 times faster



