A semi-automated AssayMAP Bravo protocol for robust Evotip loading and high-throughput analysis on the Evosep One

Evosep Biosystems, Denmark

- Agilent AssapMAP Bravo





EVUSEP

Angela McArdle, Magnus Huusfeldt, Lasse Falkenby, Jakob Bunkenborg, Laurent Rieux, Dorte B. Bekker-Jensen, Nicolai Bache

Chromatographic performance

Two protocols designed for Evotips

We developed two protocols for transferring samples from a 96 well sample plate to a box of Evotips. The first method ('full head') involved simultaneous transfer of a full set of samples from individual positions in a 96 well plate to Evotips. The second method ('single column') was designed to support user flexibility and allows for transfer of 8 samples (or a single column from a 96 well plate) at a time.



⁺Use Equilibration Holder

Layout of AssayMap Bravo deck for semi-automated loading protocols using the 'full head' or the 'single column' protocol.

Availability of AssayMAP Bravo loading protocols

Two holders are required to run the protocols; an 'Evotip' holder and an 'Equilibration' holder. The corresponding 3D print files can be found online at www.evosep.com/support/automation. The actual protocols to load on the AssayMAP Bravo can also be found here.

Excellent chromatographic performance

We extracted data for five different BSA peptides eluting across the 100 SPD method to assess the stability of retention times (RT) across 288 injections. This resulted in a median RT relative standard deviation (RSD) of 1.7 seconds. The median full peak width at half maximum (FWHM) is consistent at 3.6 seconds, which reflects the high reproducibility supported by Evotip disposable trap columns as well as the Evosep system itself.

		5	
Peptide precursor	m/z	4	Medi 3.6
TCVADESHAGCEK	488.53	4 (sec)	
LVNELTEFAK	582.32	WHN 3	
HLVDEPQNLIK	653.36	eak FV	
TVMENFVAFVDK	700.35	د 1	
LGEYGFQNALIVR	740.40	0	
			m/z 488 m/z 740 m

Chromatographic performance across the 288 injections with the 100 SPD method.

Loading protocol for Evotip Pure



Wash dry Evotips with 20 µl Solvent B and centrifuge at 800 g for 60 s. CONDITION Soak in Propanol until all the Evotips are pale white EQUILIBRATE Equilibrate soaked Evotips with 20 µl

Solvent A and centrifuge at 800 g for 60 s. Load samples on wet Evotips (20 µl in Solvent A) and centrifuge at 800 g for 60 s.

Wash Evotips with 20 µl Solvent A and centrifuge at 800 g for 60 s.

Transfer 100 µl Solvent A and centrifuge Evotips at 800 g for 10 s to keep them wet.





Robustness of quantification

Reproducible peak areas with all three methods

The reproducibility associated with manual loading was equally impressive compared to the two semi-automated protocols with coefficient of variation (CV) of 8.6% against 6.3% and 8.5% for the two Bravo protocols.



A detailed look at the two Bravo loading protocols

The reproducibility of peak area is not influenced by loading position indicated by peak area CV%. The 'full head' protocol provides slightly higher CVs than the 'single head' protocol



Reproducibility of peak areas for all detected peptides for each method (n=96).

Reproducibility of peak areas across rows and columns with both AssayMAP Bravo loading protocols.