

EVOSEP

Step-by-step guide

A complete, automated Opentrons OT-2 loading protocol

As a step towards end-to-end workflows, a simple and automation-friendly protocol for sample loading on Evtips has been developed. It makes use of a sandwich approach with defined air gaps between the liquid layers (Figure 1). This is then pushed through the Evtip with the OT-2 pipette for 100 seconds leaving the Evtip ready for injection on the Evosep One. The protocol has been converted into an easy-to-use HTML form, that generates a complete python script for use in the Opentrons app. It allows the user to load from 8 to 288 Evtips in a single run, requiring an Evtip Opentrons Loading kit (Evosep, EV1144) for each box of Evtips loaded at the time. Additionally, a solvent plate (USA Scientific, 1061-8150) and a sample plate (Eppendorf, 0030129512) is needed.

STEP-BY-STEP GUIDE

- 1 Prepare solvent plate with 18 ml solvent A and 6 ml 2-propanol for each Evtip box. The first Evtip box uses columns 1 and 12, the second Evtip box uses columns 2 and 11, and the third Evtip box uses columns 3 and 10.
- 2 Prepare sample plates, preferably minimum 30 μ l per well. 20 μ l will be loaded on the Evtip.
- 3 Prepare the remaining OT-2 deck as outlined below (Figure 1) with the Evtips placed in the appropriate adapter (Evosep, EV1144).
- 4 Go to <https://ot2-evtip-loading.evosep.com/> and select the parameters for the protocol by navigating through the options. Click 'Download'. For the Opentrons app, version 6.3.0 and higher, additional files as indicated are required to be downloaded and imported into the app.
- 5 Open the Opentrons App and go to the 'Protocols' tab. Import the generated protocol.
- 6 Open the protocol in the Opentrons app, click 'start setup', select your Opentrons and click 'Proceed to setup'. Calibrate pipettes if needed (Step 1). The protocol will now be analyzed on the robot.
- 7 Navigate to Labware setup (Step 2). Click 'Apply stored data' and run 'Labware Position Check'.
- 8 Click 'Start run' and follow the run log for current activity. Once completed, the Evtips can be transferred to the Evtip box, submerged in solvent A and stored at 4 °C until analysis.



Figure 1: Layout of OT-2 deck for automation of Evtip Pure sample loading.