

EVUSEP ENO

Software Installation Guide for CHRONOS

UM-004A

EV**U**SEP

Contents

1	Soft	ware installation for CHRONOS	. 1
	1.1	Introduction	1
	1.2	Installing Axel Semrau CHRONOS	1
	1.3	Installing the Evosep Eno CHRONOS driver	4
	1.3.1	Prerequisites:	4
	1.3.2	Installation procedure:	4

1 Software installation for CHRONOS

1.1 Introduction

This installation guide supports the installation of Evosep Eno driver when used with Axel Semrau[®] CHRONOS version 5.2.17. Evosep is not affiliated with Axel Semrau but offers an interface to enable integration.

Evosep is not responsible for the functionality, compatibility, or support of any third-party software. Integration with CHRONOS is subject to Axel Semrau's technical specifications, licensing terms, and software updates, which are outside of our control. It is the responsibility of the user or system administrator to ensure that all third-party requirements, configurations, and dependencies are met and remain supported.

Evosep assumes no liability for issues arising from changes, limitations, or failures in third-party software, nor does it provide support for software components not supplied by us. The Axel Semrau CHRONOS manual is to be consulted for full computer requirements and instructions on CHRONOS software. For use of Evosep Eno the Evosep Eno User Manual is to be consulted.

1.2 Installing Axel Semrau CHRONOS

- 1. Insert the CHRONOS USB dongle containing the CHRONOS software and license key.
- 2. Run the CHRONOS Installer application.
- 3. Click Yes on the Windows User Access Control (UAC) prompt.

EVUSEP

4. Click "Next" on the welcome dialog.



5. Accept the terms in the License Agreement and click "Next".



6. Enter the license key provided with the USB dongle and press "Next".

Axel Semrau Chronos 5.	2.17 Setup		-		\times
AXELSEMRAU	License Key To continue setup, you r	need a licens	e key.		
Please enter your lice	ense key				
Axel Semrau GmbH & Co. KG -					
	<	Back	Next >	Cano	el



7. Any desired components can be selected for install (none are required for Evosep Eno operation) and press "Next".

AXELSEMRAU	Choose Components Choose which features of Axel Se to install.	emrau Chronos 5.2.17 you war
Check the components you install. Click Next to continu	want to install and uncheck the con .e.	ponents you don't want to
Select components to insta	II: Chronos Program files Chronos Program files Example Methods ChemStation interface Xcalibur user menu Odor Control interface Purge and Trap Tasks	Description Position your mouse over a component to see its description.
	Miscellaneous Tasks	

8. Name the start menu folder destination and click "Install"

	Choose Start Menu Folde	r			
G AXELSEMRAU	Choose a Start Menu folder shortcuts.	for the Axel Semr	au Chro	onos 5.2	. 17
Select the Start Menu fold	r in which you would like to cre	ate the program's	shorta	uts. You	
an also enter a name to o	eate a new folder.				
Chronos					
Accessibility					
Accessories					Ш
Administrative Tools					
CutePDF Writer					
Dell Performance					Ц
Dropbox					
FortiClient VPN					
Ghostscript					
HP					
Java Development Kit					_
Do not create shortcut					
Semrau GmbH & Co. KG					
	< Bac	Install		Cance	-

9. Click "Finish" and the installation of Chronos has now been completed.

Axel Semrau Chronos 5.2.17 Setup -				
We sense good	Completing the Axel Chronos 5.2.17 Setu	Semrau p Wizard	I.	
chemistry!	Axel Semrau Chronos 5.2.17 has been installed on your computer.			
	Click Finish to close this wizard.			
6				
	< Back	Finish	Canc	el

EV**U**SEP

1.3 Installing the Evosep Eno CHRONOS driver

- 1.3.1 Prerequisites:
 - A. CHRONOS is installed on the computer.
 - B. The CHRONOS software is not running.
 - C. If an Evosep One Plugin is installed on the computer, the version must be newer than version 3.0.12. Otherwise, the below shown error message will appear at the start of the installation and the Evosep One driver needs to be uninstalled before the Evosep Eno driver can be installed or updated to a version later than 3.0.12 (available from evosep.com/support-zone/).



1.3.2 Installation procedure:

- 1. Connect the Evosep Eno instrument to the computer via ethernet cable. Note: Evosep Eno must be connected to the ethernet port with a static IP address.
- 2. Run the Evosep Eno CHRONOS Driver Windows installer.
- 3. Click "Next".





4. Tick the "I accept the terms in the License Agreement" checkbox and click "Next".

End-User License Age Please read the followin	reement ig license agreement carefully	EVOSEP		
Disclaimer Agreement	- Software Licer	nse		
Your use of <i>Evosep</i>	Eno Instrument Control S tributed under this license is	Software, hereafter subject to all terms and		
conditions of this lice and conditions of this	s statement, jou do not s statement, you may not use	a the software.		

5. In the Evosep+ applications window, select to install any desired applications, then click "Next".

vosep+ Select v their ins	hich Event	cations o <u>sep+</u> applications you n e in the list below.	need access to by setting	EVØSER
Evosep Chronos	+ applica s Plugin*	ation selection is also ava folder in Windows Start. Whisper Zoom	Minner Zoom m	Evosep Eno
		Will be installed on lo Entire feature will be i	cal hard drive installed on local hard drive	ive, yet inalysis of
	×	Entire feature will be u	unavailable	

6. Please read the information in "Prerequisites" window carefully, then click "Next".

😸 Evosep Eno Chronos Plugin Setup	- • •
Prerequisites Perform actions below	EVÖSEP
Please make sure the Evosep Eno instrument is powered on, connected to an Ethernet port on this computer (optionally via a network switch), and the connected network adapter is using static IP addresses!	
Back	lext Cancel



7. Click "Install" to begin the installation.



8. Click "Yes" on the Windows UAC screen, to allow the program to install the software.





- 9. During installation, the computer's ethernet configuration is checked, and one of the three cases below will occur:
 - a. If a single statically configured ethernet adapter exists, the required Evosep Eno configuration is added to that.
 - b. If multiple statically configured ethernet adapters exist, the user is presented with a selection dialog and must decide which one to use:



c. If no statically configured ethernet adapters exist, a cmd window will be displayed, stating that and to consult the manual for help:



10. When the install is completed click "Finish", to exit the installer.



Note that if you have already installed an Evosep One CHRONOS plugin, the Evosep Eno CHRONOS plugin will appear as a second instance of CHRONOS, as shown below.

