EVOSEP ENO

A Standardized Separation Tool for Proteomics

EVUSEP

Next generation insights


LC-MS based proteomics is not just about individual technologies. It's covers the entire workflow from sample preparation to data acquisition - ensuring sensitivity, reliability, and speed.

Evosep offers a standardized vet versatile solution that streamlines and optimizes the entire workflow - not just individual performance metrics. The versatility and ease-of-use benefits any type of experiment, from multiple, smaller projects to the continuous analysis of large sample cohorts.



Evosep Eno increases productivity

Evosep Eno is an easy-to-use, robust and reliable chromatographic solution for high performance and standardized proteomics, enabling seamless integration between sample preparation and LC-MS analysis in a sustainable and cost-efficient manner.




Versatility through standardized methods

Evosep Eno uses optimized methods with specified columns to ensure out-of-the-box consistent performance for simplified application support. Low overhead time ensures maximum utilization of the mass spectrometer.




Short learning curve and ease-of-use



Innovative sample introduction

Evotip Pure is a precision technology device with extremely high efficiency, integrating desalting with LC-MS sample introduction.



High efficiency trap column

Partial elution

Overload protection

No carryover

Workflow functionality


Seamless integration with end-to-end workflows Decouple sample preparation from MS availability

Innovative sample storage

Clean-up and concentration of liquid samples Loss-less sample recovery Long-term stability


Sustainable and scalable workflows

The Evotip is key for miniaturizing workflows by only having to digest just enough sample for a single injection, ensuring more cost-efficient and sustainable sample preparation.



Whisper Zoom sets a new standard

Whisper™ Zoom delivers sensitivity for highly versatile applications. Robust methods enable high-throughput analysis for inter-laboratory studies, while a short learning curve facilitates rapid method adaptation.



Robust high sensitivity workflows

The Whisper Zoom methods deliver cutting-edge performance while addressing the increasing requirement for scalability, ease-of-use and robustness.

The methods, specifically designed for standardized and high-sensitivity applications, leverage the unique properties of the Evotip Pure to sensitive storage and efficient LC-MS sample introduction.



Seamless integration

The Evosep Eno connects easily with leading mass spectrometers using optimized emitters and zero dead-volume connectors. Control is provided by integrated software drivers. **EV** SEP EN 🛡 Bruker Thermo Scientific

SCIEX

Comprehensive service coverage

Choose a tailored service plan¹ to meet your needs.

	Endurance Service Contract	Performance Service Contract	Excellence Service Contract
Troubleshooting wizards and online tools		\mathbf{V}	
Unlimited remote diagnosis and repair support	Priority	Priority	Priority
Software and method upgrades		\mathbf{V}	
Spare parts coverage	\mathbf{V}	\mathbf{V}	
Preventive maintenance (PM)	1 per year	1 per year	1 per year
Wear parts coverage		lacksquare	
Unlimited service repair visits			
Max. on-site repair response time	5 business days	3 business days	3 business days

¹ Terms and conditions apply

EV SEP

Making clinical proteomics 100 times more robust and 10 times faster

"We've got you covered, every step of the way"

Instrument

Consumables

Service

EN 🛡